

Z.O.E TYPE III TEMPORARY RESTORATIVE MATERIAL

Indications

Intermediate restorations, particularly in odontopediatrics and deciduous teeth. Protection of tooth after removal of de-mineralised tissues from cavities.

It can also be used for bases or as temporary cement, it also has the advantage of protection and healing the pulp.

Mode of Use:

Shake the bottle of powder before use.

- Take a portion of the powder in levelled scoops sufficient to produce, with the liquid, a mix of smooth and flexible consistency to fill the cavity.
- Place the powder on a slab and divide into two parts. Next divide these two parts into a further two parts each.
- On another part of the slab place 1 drop of liquid for each levelled scoop.
- With a stainless steel spatula, mix each portion of powder with the liquid thoroughly. Mixing each part for 20 - 30 seconds. Total time of mixing about 80 sec.
- Place the mixture in the previously prepared cavity; dried and with mechanical retention.
- Apply pressure to permit internal adjustment by asking the patient to bite firmly for approximately 5 minutes.
- Finally shape or adjust the entire surplus with a round bit.
- PRECAUTION: Replace the top of the bottle immediately after use to prevent evaporation and contamination, try to keep the dropper empty and in a vertical position to prevent deterioration of the bulb.

For Professional Dental use only



TECHNICAL & GENERAL Ltd. WARNINGS

Intention for use in Z.O.E. Type III Cement can be used for cementation of intermediate restorations (two to three months) of crowns and bridges and temporary not definitely inlays.

Composition:

Zinc oxide

Acrylic Polymer

Eugenol

Hardener.

The z.o.e. cement is not irritable to the skin. In case of contact with the eyes, wash immediately with abundant water for at least 15 minutes and if the irritation still persist consult a doctor. In case of contact with skin wash the area with soap and water. In case of inhalation clean the nose with water and seek medical assistance. Before the use of these products always use rubber dam.

Some patients can present allergy to the z.o.e. cement or any of its components. In this case stop its use. Do not store at temperature above 24°C. Close the packages immediately after dispensing of the product.

Keep out of reach of children.
For professional dental use only.

Contraindication: The material must not be used when the application techniques described are not possible.

Warranty: Technical & General Ltd. Warrants the product against any fault of manufacture of the material. Technical & General Ltd., does not take responsibility for the misuse of the product out of the intention for use specified or out of the instructions for use, or out of the precautions, the warnings, and indicated contraindications. Being the only obligation is the replacement of the defective product.

Information: Post-production:

If the use of these products on a patient or any one else causes some anomaly, please send detailed report to Technical & General Ltd.



CONAMCO, S.A. de C.V.

José María Bustillos No. 28 Bis, 1er y 2do piso.
Colonia Algarín, Delegación Cuauhtémoc, CP.
Ciudad de México, 06880, Mexico

EU REP

Ivonne Richaud.

Travessera de Gracia, 315, 2º-1ª 08025
Barcelona, Spain

NO REFRIGERATION NECESSARY