

CompenSate

SKU: BL LC



Light Cure Compomer Base Liner

CE 1639

Product Advantages

- Polymerizable light curing dental resin with polyacrylic acid and ionomer glass which contains zinc
- Specially formulated to compensate for composite restorations shrinkage thus resulting in a long lasting and safe restorations
- Radiopaque base lining material
- Specially formulated for use with adhesives, composites and conventional restorative materials leading to gap-free restoration.
- Chemically bonds strongly to

adhesives, composites and other resin-based materials

- Micromechanically bonds to dentin
- Prevents postoperative sensitivities and seals the dentine tubules.
- Continuous fluoride release
- Antibacterial properties due to zinc contain

Indications

- Blocking out undercuts and underfilling material for ceramic restorations, amalgam and composite.

Content

- 2 Syringes x 1.5g, 20 Tips



Sales & Information

Technical & General, 2 Albion Place, London W6 0QT United Kingdom
Tel: +44 20 8834 7575 Email: sales@technical-general.com
www.tgdent.com

Technical & General Ltd - Generated on 11.09.26