


Finger Spreaders

SKU: 0001646



Root Canal Obturation Condenser

CE 0123

Finger Spreaders

- For lateral condensation of gutta percha points
- Pointed tip
- Exceeds all ISO/ADA specifications
- High Quality Stainless Steel
- Great flexibility allowing application of high lateral pressure in curved canals
- Sizes 015-040 to match ISO standardized gutta percha points
- Color coded, ISO colors & sizes
- CC-cord handles with size and length markings
- Length of

tapered section 16mm

- Working length 25mm
- Pre-mounted silicone stoppers
- Packed in blister cards; 6 Pieces/ card in length of 25 mm
- Sterilizable up to 200C/390F

Available Sizes

Length	Size
25mm	015, 020, 025, 030, 035, 040 & assort. 015-040



Sales & Information

Technical & General, 2 Albion Place, London W6 0QT United Kingdom
Tel: +44 20 8834 7575 Email: sales@technical-general.com
www.tgdent.com

Technical & General Ltd - Generated on 11.09.26