

Flexi- Files

SKU: 1643



CE 0123

Root Canal Preparation

Flexi- Files

- Exceeds all ISO/ADA specifications
- Highly Flexible
- High Quality Special Stainless Steel Alloy
- Hand use with mounted endo stoppers
- Blister card packing 6 pieces/blister in lengths of 21 & 25 mm
- Color coded, ISO colors & sizes
- Excellent cutting ability thanks to the special Stainless Steel

Alloy with high cutting ability, precision cutting edges and unique grinding technique

- Flatter cutting angles towards the tip ensure high resistance to fracture
- Special steel alloy with Uniform structure for maximum resistance to fracture & flexibility
- The inactive tip conducts the instrument safely along the canal
- Sterilizable up to 200 C/390 F



Sales & Information

Technical & General, 2 Albion Place, London W6 0QT United Kingdom
Tel: +44 20 8834 7575 Email: sales@technical-general.com
www.tgdent.com

Technical & General Ltd - Generated on 11.09.26