

Gates Enlargers

SKU: GE



Root Enlargement of the Coronal Third

Gates Enlargers

- For enlargement of the coronal third
- Non-cutting tip with guide
- Exceeds all ISO/ADA specifications
- High Quality Stainless Steel
- Predetermined breaking point at the shank. If the instrument fractures it can be removed easily
- Ring markings
- Recommended speed 800-1200RPM
- Packed in blister cards; 6 Pieces/ card
- Sterilizable up to 200C/390F



Sales & Information

Technical & General, 2 Albion Place, London W6 0QT United Kingdom
Tel: +44 20 8834 7575 Email: sales@technical-general.com
www.tgdent.com

Technical & General Ltd - Generated on 11.09.26