


K-Reamers

SKU: 1641



CE 0123

Root Canal Preparation

K-Reamers

- Exceeds all ISO/ADA specifications
- High Quality Stainless Steel
- Hand use with endo stops
- Blister card packing 6 pieces/blister in lengths of 21, 25, 28 & 31 mm
- Color coded, ISO colors & sizes
- Excellent cutting ability thanks to Cr-Ni Steel with high cutting ability,

precision cutting edges and unique grinding technique

- Flatter cutting angles towards the tip ensure high resistance to fracture
- Special steel alloy with Uniform structure for maximum resistance to fracture & flexibility
- The inactive tip conducts the instrument safely along the canal
- Sterilizable up to 200 C /390 F



Sales & Information

Technical & General, 2 Albion Place, London W6 0QT United Kingdom
Tel: +44 20 8834 7575 Email: sales@technical-general.com
www.tgdent.com

Technical & General Ltd - Generated on 11.09.26

Available Types & Sizes

Length	Size
21mm	006, 008, 010, assort. 006-010, 015, 020, 025, 030, 035, 040, assort. 015-040, 045, 050, 055, 060, 070, 080, assort. 045-080, 090, 100, 110 & assort. 090-110
25mm	006, 008, 010, assort. 006-010, 015, 020, 025, 030, 035, 040, assort. 015-040, 045, 050, 055, 060, 070, 080, assort. 045-080, 090, 100, 110, assort. 090-110, 120, 130, 140, assort. 090-140 & assort. 120-140
28mm	006, 008, 010, assort. 006-010, 015, 020, 025, 030, 035, 040, assort. 015-040, 045, 050, 055, 060, 070, 080, assort. 045-080, 090, 100, 110 & assort. 090-110
31mm	006, 008, 010, assort. 006-010, 015, 020, 025, 030, 035, 040, assort. 015-040, 045, 050, 055, 060, 070, 080, assort. 045-080, 090, 100, 110, assort. 090-110, 120, 130, 140, assort. 090-140, assort. 120-140



Sales & Information

Technical & General, 2 Albion Place, London W6 0QT United Kingdom
Tel: +44 20 8834 7575 Email: sales@technical-general.com
www.tgdent.com

Technical & General Ltd - Generated on 11.09.26