


Paste Carriers

SKU 1645



Paste Carriers Type L

- Exceeds all ISO/ADA specifications
- Type L with safety spring
- High Quality Stainless Steel
- Metal shank
- Packed in blister cards; 4 Pieces/ card in length of 21 & 25 mm
- Color coded, ISO colors & sizes
- The cone shaped spiral of this root filler optimally utilizes the volume of the cone shaped canal thus achieving great time saving
- It prevents air bubbles inclusions due to the conical construction of

the spiral thus achieving excellent sealing

- The additional safety spiral at the end of the shank reduces the danger of breakage
- Sterilizable up to 200 C 390 F

Available Sizes

Length	Size
21mm	25, 30, 35, 40 & 25 - 40
25mm	25, 30, 35, 40 & 25 - 40



Sales and Information

Technical & General Ltd, 2 Albion Place, London W6 0QT
Tel. 020 8834 7575 Email: sales@technical-general.com
www.tgdent.com

