

## Peeso Enlargers

SKU: 1683



### Peeso Enlargers

- For enlargement of the coronal third
- Non-cutting tip with guide
- Exceeds all ISO/ADA specifications
- High Quality Stainless Steel
- Ring markings
- Recommended speed 800-1200RPM
- Packed in blister cards; 6 Pieces/ card
- Sterilizable up to 200C/390F



#### Sales & Information

Technical & General, 2 Albion Place, London W6 0QT United Kingdom  
Tel: +44 20 8834 7575 Email: [sales@technical-general.com](mailto:sales@technical-general.com)  
[www.tgdent.com](http://www.tgdent.com)

Technical & General Ltd - Generated on 11.09.26