

seBond w/Activator



Self Etching Bonding System

CE 1639

seBond w/Activator is a light-curing, self etching, self conditioning dentine/enamel Adhesive/bonding system

Indications

- Bonding system for all Light Cure or Self Cure composites, compomers and ormoceres
- For dentine and enamel. Compatible with all prevalent products

Advantages

- Latest state of the art technology (6th generation), suitable for all prevalent composites
- Very

high shear bond strength for all light or self cured restorative materials

- Easy and time saving application; Etching, priming and bonding in two steps only
- No need for etching
- Water based primer
- No post operative sensitivity
- Reliable depth of demineralization
- Long lasting restorations

Packaging

Kit1: 5ml Primer, 5 ml Bond, 5 ml Activator, 2 x 25 Micro brushes, mixing pallet



Sales & Information

Technical & General, 2 Albion Place, London W6 0QT United Kingdom
Tel: +44 20 8834 7575 Email: sales@technical-general.com
www.tgdent.com

Technical & General Ltd - Generated on 11.09.26