

tg adseal

SKU: AS01



CE 1639

Root Canal Sealer

tgadseal is an epoxy matrix resin-based root canal sealer

Indications

Permanent root canal sealer

Advantages

- Easy to mix paste-paste, automatically outputs volume units; ratio 2:1 of Base and Catalyst
- Hermetic sealing ability
- No staining to the teeth
- Insoluble in tissue fluids
- Excellent biocompatibility
- Good radiopacity

- Working time of 35mins at 23 C (73 F)
- Setting time of 45mins at 37 C (99 F)
- Does not shrink during setting
- Compatible with all gutta percha obturation techniques
- Contains no irritants and adheres to canal walls and gutta percha ensuring a permanent and excellent root canal seal

Contents

13.5g dual syringe (Base 9g, Catalyst 4.5g) , Spatula , Mixing Plate



Sales & Information

Technical & General, 2 Albion Place, London W6 0QT United Kingdom
Tel: +44 20 8834 7575 Email: sales@technical-general.com
www.tgdent.com

Technical & General Ltd - Generated on 11.09.26