

tg bluegel

SKU: EE 02

1 Syringe



4 Syringes



CE 0413

Acid etching of enamel and dentin

Indications

Acid etching of enamel and dentin

Advantages

- **tgBluegel** is an ortho-phosphoric acid etch gel which stays in place and will not run onto exposed tissues
- Dispensed from a syringe with 2.5 grams with tips, offers precise placement of gel even in relatively inaccessible areas
- **tgBluegel** washes off easily with water

Packing

EE 4K Kit of 4 syringes 8 disposable tips

EE 03 12gm syringe 10 disposable tips



Sales & Information

Technical & General, 2 Albion Place, London W6 0QT United Kingdom
Tel: +44 20 8834 7575 Email: sales@technical-general.com
www.tgdent.com

Technical & General Ltd - Generated on 11.09.26