

tg CaviLiner - LC

SKU: CV LC



CE 1639

High Quality Light Cure Cavity Liner

tgCaviliner is a high quality light cure cavity liner; hydroxide calcium phosphate in a urethane dimethacrylate resin.

Advantages:

- Excellent biocompatibility with dentine and pulp
- Light curing type cavity liner, cures extremely hard
- Insoluble in water and oral fluids
- Releases fluoride
- Radiopaque
- Easy handling

Packaging

2 x 2g Syringes, 10 disposable tips



Sales & Information

Technical & General, 2 Albion Place, London W6 0QT United Kingdom
Tel: +44 20 8834 7575 Email: sales@technical-general.com
www.tgdent.com

Technical & General Ltd - Generated on 11.09.26