

tg Cemtemp Multi

SKU: TC MP



Multipurpose Temporary Cement

Multipurpose Temporary Cement

tgCemTemp Multi- Purpose Temporary Cement is a Zinc Oxide Eugenol formulation of powder and liquid with different extracts, Balsams and Oils. When mixed, it forms a temporary cement that can be used for the purpose of:

1. Dressing/Surgical to serve as a protective covering over operated or treated areas to prevent contact with irritating or infecting substances,
2. Temporary Filling material for cavities and
3. Temporary Luting Cement for the cementation of temporary crowns and bridges.

Contents: 1 bottle 60g powder, 1 bottle 30ml liquid, 1 dropper and 1 measuring scoop



Sales & Information

Technical & General, 2 Albion Place, London W6 0QT United Kingdom
Tel: +44 20 8834 7575 Email: sales@technical-general.com
www.tgdent.com

Technical & General Ltd - Generated on 11.09.26