

tg Cemtemp NE



CE 1639

Non-Eugenol Temporary Cement

Indications

- Temporary cementation of provisional and permanent crowns bridges

Advantages

- Optimal adhesion and removal
- Does not influence or cause plasticising of resin based materials
- Easy & clean removal of excess
- Easy to mix by paste to paste type
- Thermal pulp protection
- 1.5 minutes working time after start of mixing
- 6 minutes setting time after start of mixing
- Especially suitable for patients who are sensitive to eugenol.

Contents:

Base 60g, Catalyst 20g , Mixing Pad



Sales & Information

Technical & General, 2 Albion Place, London W6 0QT United Kingdom
Tel: +44 20 8834 7575 Email: sales@technical-general.com
www.tgdent.com

Technical & General Ltd - Generated on 11.09.26