

tg ChelCream



CE 1639

EDTA Cream for Root Canal Cleaning and Preparation

Applications

- Root canal preparation with files and reamers
- Presence of calcifications in the canal
- Too narrow canals in posterior teeth
- Smear layer removal

Advantages

- Gel formula ready to use with supplied tips
- Water soluble
- Lubricates root canal and facilitates endodontic instruments penetration
- Makes shaping more efficient and easy
- High percentage

content of EDTA making cleaning of smear layer and debris even easier and efficient.

Mode of Action

- Removal of organic substances by chelation (EDTA- Ca^{2+}).
- Easier filling due to its lubricating action.
- Excellent antibacterial effect.
- Produces the cleanest canal wall after smear layer removal

Contents:

- tgChelCream 5.4ml (7.0gm) x 2, 10 disposable tips



Sales & Information

Technical & General, 2 Albion Place, London W6 0QT United Kingdom
Tel: +44 20 8834 7575 Email: sales@technical-general.com
www.tgdent.com

Technical & General Ltd - Generated on 11.09.26