

tg dam



Dental rubber dam

Product Advantages

- The main function is to isolate the tooth being treated from the bacteria and the oral cavity
- The rubber dam is held over individual teeth or groups of teeth by appropriate rubber dam clamps along the edge of the gum
- It is generally considered as non-hazardous, class 1 non-sterile, non-measuring product
- 152mm X 152mm squares of 0.18-0.23mm sheets.
- Medium thickness of green color



Sales & Information

Technical & General, 2 Albion Place, London W6 0QT United Kingdom
Tel: +44 20 8834 7575 Email: sales@technical-general.com
www.tgdent.com

Technical & General Ltd - Generated on 11.09.26