

tg LineCal Classic

SKU: LC 20



CE 0413

High Quality Calcium Hydroxide Composition

tgLineCal Classic is a high quality Calcium Hydroxide composition, rigid, Self-setting material for use in pulp capping and as a protective base-liner under dental filling materials, cements, and other base materials. It does not inhibit the polymerization of acrylic and composite restorations.

Advantages

- Self setting pulp capping agent
- Protective barrier between dentin and pulp and acid-containing cements & restorative

materials

- Suitable for use with amalgams, composites, varnishes and other permanent restorative materials
- Easily mixed and has adequate flow ability through first half of working time
- Adequately strong to resist the packing pressures
- Radiopaque
- Constant working time

Contents

tube of 13g base paste, tube of 11g catalyst paste, Applicator, Mixing pad, Red spout nozzle for base, blue Spout nozzle for catalyze



Sales & Information

Technical & General, 2 Albion Place, London W6 0QT United Kingdom
Tel: +44 20 8834 7575 Email: sales@technical-general.com
www.tgdent.com

Technical & General Ltd - Generated on 11.09.26