

tgloy Plus(+) SafeCap Capsules

SKU: LO CPLUS



Dental Amalgam Alloy Capsules

tgloy is a Premium High Silver Content (68.4%) Dental Amalgam Alloy

Among its many excellent characteristics are:

- Formulated from 99.99% pure silver, the particles or micro granules are so sized, shaped and free of contaminating surface coatings that their surfaces abrade efficiently with light triturating
- This provides a more positive and faster amalgamation than is possible with coarser grained alloys and alloys that exhibit a tendency to form surface oxides
- The scientific blend of particle sizes do not only provide ease of condensing with small-faced condensers but are also arranged to ensure an extremely dense restoration
- No fine-like dust particles, which cause high flow and variations in setting

expansion

- The clean, solid micro granules pack easily and carve with a smooth dense surface
- Strong margins can be established without breakage, as there are no coarse-jagged particles to "pull-out", leaving holes or rough surfaces
- Achieve restorations polish to the highest possible luster
- Very low static creep factor which helps prevent marginal ditching and thus increases the life and service of the restoration
- Available in both powder and capsules

Packaging

tgloy SafeCap or ScrewCap Capsules are packed in tubs of 50 Capsules each. The SafeCap Capsules come with a small plastic stick to open the capsul with.

Capsules: 1 Spill (LO 1S) , 2 Spills (LO 2S) , 3 Spills (LO 3S)



Sales & Information

Technical & General, 2 Albion Place, London W6 0QT United Kingdom
Tel: +44 20 8834 7575 Email: sales@technical-general.com
www.tgdent.com

Technical & General Ltd - Generated on 11.09.26



Committed to Excellence

Technical & General Ltd

Restoratives



Technical & General



Sales & Information

Technical & General, 2 Albion Place, London W6 0QT United Kingdom
Tel: +44 20 8834 7575 Email: sales@technical-general.com
www.tgdent.com

Technical & General Ltd - Generated on 11.09.26