

tg lube

SKU: LUB 01



All Purpose Handpiece Oil Spray

tg lube is an advanced formula universal headpiece lubricant

Advantages

- All purpose headpiece oil spray lubricant
- Specially formulated for the maintenance of all high speeds, slow speeds, angles, nosecones, head attachments and air motors
- Withstands pressure in excess of 100,000 PSI
- Withstands the high temperatures of

autoclaving to meet sterilization procedures

- Flushes and cleans debris from headpiece
- Increases appreciably turbine cartridge and bearing life
- Environmentally safe for staff and patients
- Contains trace of peppermint oil for freshness
- Universal head suitable for use with all nozzles

Packaging

500 ml Aerosol container with cap



Sales & Information

Technical & General, 2 Albion Place, London W6 0QT United Kingdom
Tel: +44 20 8834 7575 Email: sales@technical-general.com
www.tgdent.com

Technical & General Ltd - Generated on 11.09.26