

tg lutoglass



CE 0413

Glass Ionomer Cement Type I Luting

tg lutoglass is a glass ionomer Type I Luting cement specifically formulated as the ideal ultra-thin glass ionomer crown and bridge luting cement

Indications

For cementing crowns, bridges & posts

Advantages

- Specifically formulated to attain high adhesion to tooth &

restoration

- Strong bond to dentine and fluoride ion release assures a long lasting bonding and strength integrity
- High diametric tension to withstand masticator forces
- Universal tooth shade
- Ultra-thin film thickness: 15 microns. Radiolucent

Contents

20 gm Powder, 15ml Liquid and measuring scoop



Sales & Information

Technical & General, 2 Albion Place, London W6 0QT United Kingdom
Tel: +44 20 8834 7575 Email: sales@technical-general.com
www.tgdent.com

Technical & General Ltd - Generated on 11.09.26