

tgmicrohybrid Composite

SKU MH



tgmicrohybrid (Single Syringe)
Micro Hybrid Light Cure Composite

Indications

- Universal Composite for restorations of all cavity classes
- Indirect inlays, onlays and veneers

Advantages

- Direct posterior and anterior restorations of all cavity classes I-V
- Excellent polishing and sculpting properties
- Optimal mechanical properties, maximum wear resistance and fracture strength
- Low polymerization

shrinkage and excellent aesthetics

Shades

- Available shades to choose from are: A1, A2, A3, A3.5, A4, B2, B3, C1, C2, C3, D2, D3. Kits are available upon request.

Options

- The light cure composite syringes can be purchased separately or in kits. The kits are in 7, 5 and 2 syringes. The shades in the kits can be chosen to suit your requirements

Contents:

1 Syringe: 1 x 4.5 gm Syringe in smart carton box



Sales and Information

Technical & General Ltd, 2 Albion Place, London W6 0QT
Tel. 020 8834 7575 Email: sales@technical-general.com
www.tgdent.com

