

tg orthoglass

SKU: OD



CE 0413

Glass Ionomer Type I Orthodontic

tgOrthoglass is a Type I Orthodontic Glassionomer Cement

Indications

For cementation use in Orthodontics

Advantages

- Strong bond to metal, band and enamel
- Gives sufficient working time for multiple placements
- Universal tooth shade
- Ultra-thin film thickness: 15 microns. Radiolucent

Contents:

20 gm Powder 15ml Liquid and measuring scoop



Sales & Information

Technical & General, 2 Albion Place, London W6 0QT United Kingdom
Tel: +44 20 8834 7575 Email: sales@technical-general.com
www.tgdent.com

Technical & General Ltd - Generated on 11.09.26