

**tg paste**

SKU: TP



CE 1639

### Calcium Hydroxide with Barium Sulfate

#### Calcium Hydroxide with Barium Sulfate

**tgpaste** is a water soluble temporary root canal filling material consisting of Calcium Hydroxide with Barium Sulfate

#### Indications

- Pulp capping and pulpotomy
- Leaking canal
- Formation of a hard tissue barrier (Apexification)
- Temporary root canal filling
- Internal & external root resorption

#### Advantages

- Easy cleaning and removing with good

#### water solubility

- Excellent antibacterial effect and radiopacity
- Premixed paste
- Syringe type for easy delivery of the paste into the canal
- Disposable tips for prevention of cross-contamination

#### Contents

- 1 Syringe Kit: 2.2g paste, 20 Disposable tips, one ring rotator for direction control of the tip
- 2 Syringes Kit: 2.2g paste, 20 Disposable tips, one ring rotator for direction control of the tip



#### Sales & Information

Technical & General, 2 Albion Place, London W6 0QT United Kingdom  
Tel: +44 20 8834 7575 Email: sales@technical-general.com  
www.tgdent.com

Technical & General Ltd - Generated on 25.07.26