

tg Pit & Fissure LC

SKU: PF LC



CE 0413

Visible Light Curing Sealant Opaque

tgPit & Fissure LC

Pit & Fissure Visible Light Curing Sealant Opaque
tgPit & fissure LC sealant is a light curing bi-Gma based sealant resin formulated specifically for the sealing of small fissures, pits and occlusal surfaces for the purpose of caries prophylaxis. The flowable sealants have thixotropic effects making them excellent for penetration into the tight adaptations with high bond strength to the enamel achieving high wearing resistance.

Advantages

- Classical light cure sealant
- Opaque with adjusted viscosity
- High bond strength
- Highly resistant to wear
- Excellent penetration characteristics
- Shelf life of 3 years minimum

Additional Instructions

- Shake sealant syringe very well before each use

Contents:

Light Cure Kit: 3 x 2g Sealant syringes, 2.5g Etchant Syringe, 5 Microbrushes & 8 tips



Sales & Information

Technical & General, 2 Albion Place, London W6 0QT United Kingdom
Tel: +44 20 8834 7575 Email: sales@technical-general.com
www.tgdent.com

Technical & General Ltd - Generated on 11.09.26