

tgPit & Fissure

SKU PF SC



Pit & Fissure Self Curing Sealant Tinted(Clear/Opaque).

tgPit & fissure SC sealant is a self curing bi-Gma based sealant resin formulated specifically for the sealing of small fissures, pits and occlusal surfaces for the purpose of caries prophylaxis. The flowable sealants have thixotropic effects making them excellent for penetration into the tight adaptations with high bond strength to the enamel achieving high wearing resistance.

Advantages

- Classical self cure sealant
- Tinted with adjusted viscosity
- High bond strength
- Highly resistant to wear
- Excellent penetration characteristics
- Shelf life of 3 years minimum

Additional Instructions

- Shake bottle before each use

Contents:

Self Cure Kit: 3ml Sealant Base, 3ml Sealant Catalyst, 2.5g Etchant syringe, 5 Mixing Wells, Mixing pad, Brush Holder & 20 Mini Brushes



Sales and Information

Technical & General Ltd, 2 Albion Place, London W6 0QT
Tel. 020 8834 7575 Email: sales@technical-general.com
www.tgdent.com

