

tg putty Handmix

SKU: PT



highly proven Impression Material based on VPS

Suitable for first impressions and other things.

Available in **standard, fast and soft versions.**

Indications:

Crown and Bridge impressions, Inlay and Onlay impressions

Advantages:

- Easy dosage and easy removal from mouth
- Excellent dimensional stability and excellent recovery after deformation
- Outstanding kneadability
- Exceptional when used for two phase impressions with tgPutty light/fast

Packaging: 2 x Jars of 262ml Base and Catalyst paste with 2 spoons



Sales & Information

Technical & General, 2 Albion Place, London W6 0QT United Kingdom
Tel: +44 20 8834 7575 Email: sales@technical-general.com
www.tgdent.com

Technical & General Ltd - Generated on 11.09.26