

tg Resincem

SKU: DC 01

1 Syringe

Automix



Dual Cure Resin Cement

4 Syringes



Dual Cure Resin Cement

CE 0413

tgResincem is a Dual Cure Resin Cement for permanent cementation of valuable restorations such as composites, crowns, bridges, veneers, inlays and onlays.

Advantages

- Combines the merits of glassionomers and composites
- Ideal for the cementation of Maryland bridges, veneers and facings and full coverage restorations
- Universal use and suitable for ceramic works

- Its formula blends perfectly with any restoration
- Flowable but not slumpy thus providing low film thickness of less than 15 µm
- The radiopaque barium glass filler with the resin matrix guarantees outstanding abrasion resistance
- High mechanical strenght
- Continuous fluoride release

Contents

Kit of: 2 x 2.3g syringes base, 2 x 2.3g syringes catalyst, 8 tips



Sales & Information

Technical & General, 2 Albion Place, London W6 0QT United Kingdom
Tel: +44 20 8834 7575 Email: sales@technical-general.com
www.tgdent.com

Technical & General Ltd - Generated on 11.09.26