

tg restoglass

SKU: GI RES 21



CE 0413

Glass Ionomer Cement Restorative Type II

tgrestoglass is a Glass Ionomer restorative Type II new improved version of tg widely acclaimed Glass Ionomer Cement Type II. It is ideal for glass ionomer restorations.

Indications

For class III and V restorations
tunnel preparations, root caries, and deciduous teeth preparations and underfillings.

Advantages

Excellent aesthetics: opacity comparable with composites
Pronounced condensable setting phase
High abrasion resistance
Very low early solubility
Very good bonding to dentin and enamel
Universal tooth shade; Shades 21 & 22

Contents: 20 gm Powder, 15ml Liquid, 10g Varnish and measuring scoop



Sales & Information

Technical & General, 2 Albion Place, London W6 0QT United Kingdom
Tel: +44 20 8834 7575 Email: sales@technical-general.com
www.tgdent.com

Technical & General Ltd - Generated on 11.09.26