

tg Restotemp

SKU: ZO T3



CE 0413

Temporary Restorative Material

tgrestotemp Temporary Restorative Material is a Z.O.E. Type III Intermediary Restorative Material

- For Odontological intermediate restorations, particularly in odontopediatrics
- For use with rampant carries where the de-mineralized tissue should all be removed from the cavity as fast as possible to reduce the concentration of cariogenic bacteria and prevent the processes

that produces carries. During this time the teeth ought to be protected by tgrestotemp which is a restoration of intermediate duration

- **tgrestotemp** can also be used for bases or as temporary cement. It also has the advantage of protection and healing the pulp

Contents:

(A) 38gm Powder (B) 14ml Liquid (C) Dropper and measuring scoop



Sales & Information

Technical & General, 2 Albion Place, London W6 0QT United Kingdom
Tel: +44 20 8834 7575 Email: sales@technical-general.com
www.tgdent.com

Technical & General Ltd - Generated on 11.09.26