

tg sealer

SKU: RCS 01



CE 0413

ZOE Root Canal Sealer

tgsealer is a Zinc Eugenol Root Canal Sealer

Indications

ZOE Root Canals sealer

Advantages

- Does not irritate periapical tissues
- Anti-inflammatory and contains antiseptic and antichloristic additives for the relief of post operative pains following root canal procedure
- Does not decolorize the structure of

the tooth

- Due to its content of Barium Sulphate, it is radiopaque and gives an excellent contrast for better diagnostic control
- It's high adhesion properties seals hermetically
- Fine flowing consistency for easy application
- Long workable properties

Packaging

20ml Zinc oxide Powder, 15 ml Eugenol Liquid, Dropper and a measuring scoop



Sales & Information

Technical & General, 2 Albion Place, London W6 0QT United Kingdom
Tel: +44 20 8834 7575 Email: sales@technical-general.com
www.tgdent.com

Technical & General Ltd - Generated on 11.09.26