

tg tempolux

SKU: LX



CE 1639

Automix Temporary Crown & Bridge Material

Indications:

Temporary crowns and bridges, temporary inlays and onlays, long term temporaries and temporary veneers.

Advantages:

- exceptional flexural strength, highest pressure and abrasion resistance;
- exactly fitted and unstressed due to the slightest polymerisation shrinkage;
- minimal polymerisation temperature during setting;
- defined elastic phase for reliable removal of the provisional;
- brilliant, translucent colours for an exceptional natural aesthetic and fluorescence even in different lighting

conditions;

- absolutely colour stable for complete time of wearing, biocompatible, tasteless;
- colour-stable and publishable to a high shine, providing for a natural aesthetic. Its colouration can be individually adjusted, compatible with UV-light-polymerising composites;
- low polymerisation temperatures prevent pulpal and gingival irritation and damage;
- for direct application in Automix2 cartridge system, ideal mixing ratio of 10 : 1;

Contents:

Cartridge 50ml/76gm, 10 mixing tips blue. Available in shades A 1, A2, A3, A3.5 and B1



Sales & Information

Technical & General, 2 Albion Place, London W6 0QT United Kingdom
Tel: +44 20 8834 7575 Email: sales@technical-general.com
www.tgdent.com

Technical & General Ltd - Generated on 24.07.26