

tgzoe I

SKU: ZO T1



CE 0413

ZOE Temporary Cement Type I

ZOE Temporary Cement Type I

tgzoe I is an Improved Zinc Oxide - Eugenol Cement for temporary Cementation and filling

Indications

Recommended for temporary fillings and cementation

Advantages

- May be used as a sedative base under metallic restorations
- High tensile and compressive strength
- The setting time range is from 6 to 8 minutes
- Protects the pulp from thermal shock and possible irritation from freshly mixed cements

Contents: 20 gm Powder 15ml Liquid and measuring scoop



Sales & Information

Technical & General, 2 Albion Place, London W6 0QT United Kingdom
Tel: +44 20 8834 7575 Email: sales@technical-general.com
www.tgdent.com

Technical & General Ltd - Generated on 11.09.26