Committed to Excellence
Technical & General Ltd

Cemtemp NE



CE 1639

Non-Eugenol Temporary Cement

Indications

• Temporary cementation of provisional and permanent crowns bridges

Advantages

• Optimal adhesion and removal
• Does not influence or cause plasticising of resin based materials
• Easy & clean removal of excess
• Easy to mix by paste to paste type
• Thermal pulp protection
• 1.5 minutes working time after start of mixing
• 6 minutes setting time after start of mixing
• Especially suitable for patients who are sensitive to eugenol.

Contents:

Base 60g, Catalyst 20g , Mixing Pad

Sales & Information
Technical & General, 2 Albion Place, London W6 0QT United Kingdom
Tel:+44 20 8834 7575 Email: sales@technical-general.com
www.tgdent.com