Committed to Excellence
Technical & General Ltd

# tgChelCream




CE 1639

EDTA Cream for Root Canal Cleaning and Preparation

**Applications**

- Root canal preparation with files and reamers
- Presence of calcifications in the canal
- Too narrow canals in posterior teeth
- Smear layer removal

**Advantages**

- Gel formula ready to use with supplied tips
- Water soluble
- Lubricates root canal and facilitates endodontic instruments penetration
- Makes shaping more efficient and easy
- High percentage content of EDTA making cleaning of smear layer and debris even easier and efficient.

**Mode of Action**

- Removal of organic substances by chelation (EDTA-$Ca^2$).
- Easier filling due to its lubricating action.
- Excellent antibacterial effect.
- Produces the cleanest canal wall after smear layer removal

**Contents:**

- **tg**ChelCream 5.4ml (7.0gm) x 2, 10 disposable tips

**Sales & Information**
Technical & General, 2 Albion Place, London W6 0QT United Kingdom
Tel:+44 20 8834 7575 Email: sales@technical-general.com
www.tgdent.com