Committed to Excellence
Technical & General Ltd

tg implacem

SKU: IC SP



CE 1639

Suitable with all Material Combinations

tgimplaCem is the most up-to-date developed luting cement, dual curing, semi-permanent, resin based occlusally stable specially developed for implant-retained restorations and all super structures such as crowns and bridges. tgimplaCem is suitable with all material combinations.

Advantages

• Passes through a plastic phase before final setting thus allowing easy, convenient &clean removal of excess cement.

• Very low displacement resistance and thus extremely low film thickness.

• Due to highly cross-linked, ductile polymer structure, high compressive strength andhermetic sealing are achieved

• Luted restorations are effortlessly removed without damage with standard devices(Crown Butler, Corona Flex).

• Excess material can easily be removed in large pieces (no crumbling) during plasticphase.

• Compatible with all conventional material combination, i.e.Zirconium and metalframeworks, ceramics, titanium.

• Antibacterial, eugenol-free, odourless and tasteless; natural opaque colour for screeningmetallic abutments.

Packing & Contents: 1x5 ml Syringe, 10 Tips

Sales & Information
Technical & General, 2 Albion Place, London W6 0QT United Kingdom
Tel:+44 20 8834 7575 Email: sales@technical-general.com
www.tgdent.com